

Release Notes **Tobii Dynavox Communicator 5**

Introduction

Communicator 5.0.1 is a minor update to Communicator 5 that fixes 3 problems facing Communicator 5; Facebook Mandated API migration, Problems saving TC4 page sets, and weaknesses related to sending/receiving text messages. This Update is compatible with all I-Series devices, Windows Tablets, and PCs.

System Requirements

- Windows 7, Windows 8, Windows 8.1 and Windows 10
- 2.0GB RAM
- 3.0GB Hard Drive
- 1.8 Ghz processor

Customer support: For more information or support, please visit www.TobiiDynavox.com.

Release Version 5.0.1

Release date: 2015-09-22

Improvements

- Improved reliability when sending or receiving text messages through a Bluetooth phone (CFIVE-6573)
- Updated Facebook to use the new API 2.0 (CVI-17)
- Occasionally unable to Save Communicator 4 Page Set in Communicator 5 (TDCOM-258)

Support for Your Tobii Dynavox Device

Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: www.TobiiDynavox.com or www.myTobiiDynavox.com.

Contact Your Sales Representative or Reseller

For questions or problems with your product, contact your Tobii Dynavox sales representative or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit www.TobiiDynavox.com/contact