

Release Notes **Tobii Dynavox Communicator 5**

Introduction

Tobii Dynavox Communicator 5 is a software package designed to help individuals with communication disabilities communicate more efficiently. Communicator 5 converts text and symbols into clear speech, and offers easy-to-use tools for computer access, e-mail, text messaging, telephony, environmental control and more.

System Requirements

- Windows 7, Windows 8, Windows 8.1 and Windows 10
- 2.0GB RAM
- 3.0GB Hard Drive
- 1.3 Ghz processor

Customer support: For more information or support, please visit www.TobiiDynavox.com.

Release Version 5.1

Release date: 2015-12-09

New Features

- The new Abbreviation Expansion feature allows users to access their phrases faster than ever.
- There are new Alarm and Go To Desktop page sets
- The user can set a Different voice to be used for speech output and feedback sounds
- The Combine 2 IR commands action has been updated to allow the user to also Combine 3 IR commands
- Added Keyboard support for Croatian, Czech, Hungarian, Polish, Slovak, Slovenian, & Ukrainian
- Communicator 5 is now available in Finnish, Brazilian Portuguese, & Japanese

Improvements

- Communicator leaks memory when loading cddx files (TDCOM-327)
- Links refer to www.pagesetcentral.com (TDCOM-311)
- Importing TC4 E-Mail into C5 may cause a crash on startup (TDCOM-256)
- Unexpected results when injected page has the same name as a page in the main page set (TDCOM-193)
- Skype OSK Page Set cannot place calls (TDCOM-294)

Release Version 5.0.1

Release date: 2015-09-22

Improvements

- Improved reliability when sending or receiving text messages through a Bluetooth phone (CFIVE-6573)
- Updated Facebook to use the new API 2.0 (CVI-17)
- Occasionally unable to Save Communicator 4 Page Set in Communicator 5 (TDCOM-258)

Support for Your Tobii Dynavox Device

Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: www.TobiiDynavox.com or www.myTobiiDynavox.com.

Contact Your Sales Representative or Reseller

For questions or problems with your product, contact your Tobii Dynavox sales representative or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit www.TobiiDynavox.com/contact