

Release Notes

DynaVox Compass Software, Version 2.1.5

This document summarizes the new features and enhancements and corrected issues associated with the release of version 2.1.5 of the DynaVox Compass Software.

New Features and Enhancements for 2.1.5

- Compass is now fully integrated with the Tobii I-12 and I-15 devices. Specifically, Compass can now leverage the full support of eye tracking and IR-based environmental control.
- Full support for Spanish has been introduced with this release, including:
 - A Spanish language choice in the Interface Settings menu,
 - Spanish strings throughout the system menus,
 - Spanish word prediction and word morphology,
 - Spanish TTS voices, and
 - Spanish versions of the All Access and Navigator page sets.
- IR learning is introduced for the T10
- The myDynaVox account creation process is now fully supported within the Compass software enabling new users to be up and running without ever having to open a web browser.
- Seamless interaction with myDynaVox allows the Compass software to present various menus such as importing pages, restoring backups, and browsing voices for download without requiring that a user logs in to their own myDynaVox account.
- Upgrading Compass no longer requires a myDynaVox account.
- Several changes have been made to the Compass software to comply with the new PDAC requirements. This ensures Medicare funding compliance for Compass on the I-12 and I-15.
- A significant bug that prevented upgrading the 2.0 Compass software has been addressed.
- Bug fixes and stability improvements to the Compass software.
- Introduced an improved and consistent user interface throughout all of the Compass access method settings menus.

New Features and Enhancements for myDynaVox

- The myDynaVox single account restriction has been removed. This enables a user to attach Compass to another myDynaVox account without deactivating (detaching) from the first account.