

# Gateway Version 1.3 Release Notes

## Added Dictionaries to the Toolbar for Child 20, Adult 20, 30, 40, and 60

- To enable quick access to fringe vocabulary, each new Gateway Dictionary contains all of the Folder Categories for the specified pageset. A user may select the Dictionary that corresponds with their pageset or a more advanced Dictionary to gain access to additional fringe vocabulary.

## Improved Scanning in all Pagesets

- All pages and Tabs have been optimized for scanning. The Open Pop-up Exclusively function is used where indicated to ensure that only the targeted page is scanned.
- Symbol alignment has been improved on pages to offer the correct scan pattern.
- Group scanning has been added where developmentally appropriate to enhance the rate of access to a targeted location.

## Common Core

- A Common Core page has been added to Gateway 30

## Phrase Prediction

- Phrase Prediction has been added to the Phrases Tab in the Phrases Tabset. Affects Child Gateway Pro and Gateway Pro.
- Added Phrase Prediction to ScanPro, including two phrase prediction buttons located beneath the word prediction buttons. Affects both the child and adult versions of ScanPro.

## Character Prediction

- Added Character Prediction to the child and adult versions of ScanPro.

## TV Controls

- Basic TV control functions are now available in the ECU Folder on the Time to Relax Tab. Affects Child 20 and Adult 20.
- Some TV controls have been added to the Watch TV folder on the time to Relax tab. Affects Gateway 30.
- TV/Remote folder has been added to the Time to Play tab. Affects Child 40.
- Added TV remote operations to the ECU/Remote folder on the Time to Relax Tab. Affects Teen/Adult 40 and Text 40.
- TV remote functions added to TV folder on the Time to Relax Tab. Affects Child 60, Teen/Adult 60 and Text 60.

## Corrections

- Added links within the Things folder. Affects child and adult versions of ScanPro.
- Added past tense for the Drive Semantic Power Strip. Affects Gateway 40 Teen/Adult, Child 60, Teen/Adult 60 and 60 Text.
- Corrected several typographical errors for the words - fruit, avocado, and cantaloupe.

## Support for Your Tobii Dynavox Device

### Get Help Online

Before contacting us, see the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support page online at: [www.tobiidynavox.com](http://www.tobiidynavox.com).

### Contact Your Reseller

For questions or problems with your product, contact Customer Support at your Tobii Dynavox Reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For our Resellers contact details, see [www.tobiidynavox.com](http://www.tobiidynavox.com)