

Release Notes **Tobii Dynavox Snap + Core First** for the **Windows version**

Release Version 1.1

Release date: 2017-10-09

New Features

- Import, export, share, and Sync your content.
- Integration with myTobiiDynavox for import, export, share, and Sync.
- Support for Windows Control 2, EyeMobile Plus, and PCEye Plus.
- New Message Window options in Page Set > Style:
 - font size
 - text color
- Added “Delay between selections” setting for Touch Exit and Touch Enter.
- Added the action “Launch Communicator.”
- Updates to Visit functionality, including three new navigation actions:
 - Cancel Visit
 - Extend Visit
 - Return from Visit
- Updated functionality for the free version of the app:
 - Companion app for Tobii Dynavox device users
 - SLPs can obtain the full functionality of the app with their valid ASHA code
 - Those not covered by the instances above may obtain a 30-day free trial of the speaking app (via myTobiiDynavox)
 - In-app purchase of speech
- Choice of content language/culture through the Setup Wizard.
 - New content for the U.K.

Improvements

- Updated Setup Wizard interface.
- Performance improvements.
- Bug fixes.

Support for Your Tobii Dynavox Device

Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: www.TobiiDynavox.com or www.myTobiiDynavox.com.

Contact Your Sales Representative or Reseller

For questions or problems with your product, contact your Tobii Dynavox sales representative or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit www.TobiiDynavox.com/contact