



Tobii Dynavox Snap + Core First Content Release Notes

Release Version 1.1

Release date: 2017-10-09

New Features

- Added the Bathroom topic and all associated pages.

Improvements

- 800+ symbol changes throughout the Page Set makes symbols more consistent across common vocabulary.
- Common QuickFires words and simple expressions have been added to the appropriate word lists, making those words and phrases accessible for writing.
- All toolbar buttons (except the Back and Word Forms buttons) now have the Cancel Visit action. This prevents unexpected visit page returns.
- Added "am" and "are" to the Actions word list.
- Added missing "u" key to the 6x6 navigation button grid size for the grid matching keyboard.
- Full-sized keyboards are now 8x10 for keyguard compatibility.
- Reordered the following word lists:
 - Clothing — clothing items appear first in alpha order, followed by the colors in alpha order.
 - Money — coins and dollars appear in ascending amounts, followed by the remaining words in alpha order.
 - Seasons — weather words appear first in alpha order, followed by the seasons in chronological order.
 - Minor ordering fixes to multiple word lists.
- Added missing symbols and removed extra blank buttons in several word lists.