

Snap Scene Release Notes

Release Version 1.2

Release date: 2016-9-15

Enhancements and Improvements to Snap Scene

- Snap Scene is now available in simplified Chinese.
 - To install the software in a language other than English, the device operating system must be set to the desired language at the time of installation.
 - Changing the operating system language after Snap Scene is installed will localize the labels in Snap Scene, but not the default content.
- Improved UI to allow editing and playing scenes to happen more naturally.
- Improved hotspot label animation.
- Added Show Hotspot Labels setting option “On” to show hotspot labels at all times.
- A High-contrast Hotspot Label setting is now available.
- The Show Hotspots setting now has the “Fade Out” option, which shows the hotspots when a scene is first opened, then fades them away.
- Bug fixes.

Release Version 1.1

Release date: 2016-6-15

Enhancements and Improvements to Snap Scene

- Snap Scene is now available in the following languages: English, French, German, Danish, Dutch, Norwegian, Spanish, and Swedish.
 - To install the software in a language other than English, the device operating system must be set to the desired language at the time of installation.
 - Changing the operating system language after Snap Scene is installed will localize the labels in Snap Scene, but not the default content.
- The Windows platform is now supported, including open I-Series devices.
- Snap Scene is now available on open Tobii Dynavox T-Series tablets.
- On T-series devices, scanners may now use the App Switcher functionality to change between Snap Scene and Compass v2.5.
- Bug fixes.

Support for Your Tobii Dynavox Device

Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: www.TobiiDynavox.com or www.myTobiiDynavox.com.

Contact Your Sales Representative or Reseller

For questions or problems with your product, contact your Tobii Dynavox sales representative or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit www.TobiiDynavox.com/contact