

## Release Notes **Tobii Dynavox WordPower**

### Release Version 3.0

**Release date:** 2016-11-30

#### Improvements

- Improved alignment with Tobii Dynavox keyguards
- Miscellaneous vocabulary updates and improvements

#### New Features

- Addition of Scanning WordPower vocabularies:
  - Scanning WordPower-160
  - Scanning WordPower-100
  - Scanning Picture WordPower-Phrases and Categories
  - (WordPower-20 Touch & Scan)
- Addition of Remote Control (Environmental Control features)
- Addition of Tobii Dynavox Games
- Addition of Tobii Dynavox Calculator
- Added buttons that link to popular Android applications (including text messaging option through the “Send Text” app with an Android phone)

### Support for Your Tobii Dynavox Device

#### Get Help Online

See the product-specific Support page for your Tobii Dynavox device. It contains up-to-date information about issues and tips & tricks related to the product. Find our Support pages online at: [www.TobiiDynavox.com](http://www.TobiiDynavox.com) or [www.myTobiiDynavox.com](http://www.myTobiiDynavox.com).

#### Contact Your Sales Representative or Reseller

For questions or problems with your product, contact your Tobii Dynavox sales representative or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training. For contact details, visit [www.TobiiDynavox.com/contact](http://www.TobiiDynavox.com/contact)